

Art And Science

ELEGANCE & DISTINCTION Westar introduces affordable luxury. Always famous for reliability and efficiency. Westar adds beauty and comfort to the dental operator.

❖ 5000 Ultra Stools



“Elegance and Distinction - Westar introduces affordable luxury. Always famous for reliability and efficiency, Westar adds beauty and comfort to the dental operator. The classic look of contoured lines and leather will move you to the forefront of modern operator design.”

WESTAR ❖

❖ 2000 Economy Stools



Westar Medical Products, Inc.

Phone: 360/435-8900 ❖ Watts: 800/499-3945 ❖ Fax: 425/328-1544

18930 59th Ave. N.E. ❖ Arlington, WA 98223

Email: westar_medical@msn.com ❖ Website: www.westarmed.com





D5000



A5000



D2000E



A2000E



5000 Ultra Stools

2000-104-1	D5000 Ultra Doctor Stool	\$652
2000-105-1	A5000 Ultra Assistant Stool	\$872

Seating is adjustable to multiple positions with an easy to reach paddle, reducing fatigue

Lumbar support for both the doctor and assistant stools are large and adjust vertically as well as horizontally

Doctor's seat is large and comfortably contoured reducing pressure points

Assistant torso support micro adjusts vertically and ratchets easily into position

The body support has special high density foam for added comfort

Five point 21" base - powder coated cast aluminum

Dual wheel, non-marring casters

2000 Economy Stools

2000-103-1	D2000E Economy Doctor Stool	\$404
2000-102-1	A2000E Economy Assistant Stool	\$674

Doctor stool back adjusts vertically and horizontally

Assistant torso support adjusts vertically and swings into position

Height adjusts easily with the touch of a paddle

Assistant stool has an adjustable foot ring

Five point 22" plastic base

Dual wheels barrel casters

Options and Parts

"D" Seat in lieu of round	8003-199	\$152
Ultraleather upholstery upgrade	1006	\$186
Assistant Cylinder (Ultra)	8000-017-1	\$96 or in lieu of
Doctor Cylinder (Ultra)	8000-018-1	\$96 or in lieu of
Assistant Cylinder (Economy)	8000-017	\$96 or in lieu of
Doctor Cylinder (Economy)	8000-018	\$96 or in lieu of
Adjustable foot ring	8002-214	\$108
Casters for plastic base (specify color)	8000-002	\$10
Casters for metal base - non-marring	8000-033	\$12

(Economy stools have black bases)