


ARTISTRY IN ENGINEERING



Constellation with Assistant Delivery System

**Westar Medical Products, Inc.**


Phone: 360/435-8900 ❖ Watts: 800/499-3945 ❖ Fax:  
425/328-1544 18930 59th Ave. N.E. ❖ Arlington, WA 98223  
Email: [westar\\_medical@msn.com](mailto:westar_medical@msn.com) ❖ Website: [www.westarmed.com](http://www.westarmed.com)

*Nothing in this world is so powerful as an idea whose time has come.*

--- Victor Hugo

**Constellation with Contoured Cuspidor Delivery System**

6000 \$5356



*(Shown with optional features: Delivery mounted Touchpad and WeStarlight VIII)*

*Shown on cover:*


**Constellation Delivery System with Assistant Package**

(Sold separately)

2001 Plush with Trim Patient chair	2000-401T	\$5848
D5000 Ultra Doctor Stool	2000-104-1	\$ 652
A5000 Ultra Assistant Stool	2000-105-1	\$ 872
Ultraleather upcharge - Stool	1006	\$ 186 ea
WeStarlight VIII	2000-245	\$ 964

**Constellation with Assistant Package Delivery System**

6010 \$4730



*(Shown with optional features: Delivery mounted Touchpad and WeStarlight VIII)*

**International Package**

2000-020 Cuspidor \$7920

*Not shown:*

2000-019 Assistant \$7520



**Professional Cuspidor Delivery System**


6000PC Cuspidor \$2984

*Not shown:*

6000PA Assistant \$2700

Westar Chair Adapter 2000-226 \$ 284

*(Needed if ordering Westar Chair)*



(Sold separately on Professional deliveries)

WeStarlight VIII	2000-245	\$ 964
45" Light Pole (n/c if in lieu of 27" pole)	2000-214	\$ 124
1.375" Pole Adapter	4002-071	\$ 28

**International Package**

Package includes:

- 2001E Economy Patient Chair
- D2000E Economy Doctor Stool
- Two-handpiece automatic Control
- Pole mounted arm and flex without air brake
- WeStarlight VIII with High and Low intensity settings
- One each 1/4" Air/water quick disconnects
- Pole mounted Cuspidor with Saliva Ejector vacuum (2000-020)

Or

- Telescoping arm with Saliva Ejector vacuum (2000-019)

(Sold as Complete package only)

*(Shown with optional features: Delivery mounted Touchpad, Stainless Instrument Tray, and WeStarlight VIII)*

**Standard Features:**

**Delivery Systems:**

- Three handpiece automatic valves  
*(Int'l - two handpieces)*
- Tri-Block valve *(Duo-block on 2 handpiece systems)*
- Master on/off
- Water coolant adjustment
- Three way syringe with autoclavable tip
- Non-retracting water coolant
- Handpiece purge toggle
- Asepsis Midwest 4 hole tubing
- Standard foot control
- Utility Center
- Regulators, Filters, and Air/Water shut off valves
- Electronics cable pre-wiring (2 conductor)

**Features included with**

**Constellation Delivery systems:**

- Flex arm with pneumatic brake
- Contoured Ridged Arm
- Westar chair adaptor for 2" pole
- 45" Light pole with light adapter bushing
- Stainless instrument tray

**Features included with Excursion Delivery systems:**

- Front Mount Left/Right rotation
- Clean water system *(bottle water only)*

**Features included with Professional and**

**International Delivery Systems:**

- Flex arm without air brake
- Pole Clamp Ridged Arm
- 27" Pole-Professional, 45" Pole-International
- Westar chair adapter *(included with International, purchase separately for Professional)*

**Features included with**

**Constellation Cuspidor Systems:**

- Contoured porcelain bowl
- HVE and SE with autoclavable valve
- Solids collector
- Bowl flush and cup fill - Timed
- Water quick disconnect
- Clean water system *(bottle water only)*
- Water Master Shut-off valve

**Features included with Professional and**


**International Cuspidor Systems:**

- Round porcelain bowl
- HVE and SE with autoclavable valve
- Solids collector
- Bowl flush and cup fill - Manual
- Water Master Shut-off valve

**Features included with Assistant Systems:**

- Telescoping assistant arm
- HVE and SE
- Solids collector
- Water quick disconnect *(Optional feature on Professional or Excursion)*
- Clean water system *(bottle water only) (Optional feature on Professional)*

**Excursion Delivery System**



**Euro Head**

Add \$580 to delivery price



	<u>Part #</u>	
Excursion Delivery with Nautica Control	5030	\$5110
Excursion Delivery with 6000 Esense Control	6030	\$4832
Excursion Delivery with Euro Control	5037	\$5412
Excursion Light and Post (purchased with unit)	5032-245	\$1496
Excursion Light Post only (purchased with unit)	5032	\$532
Excursion Light Post (not purchased with unit-includes Front Mount)	5032-169	\$836
Excursion Monitor Pole (purchased with unit)	5034	\$852
Front Mount for Poles purchased without unit	4000-169	\$306
Flat Screen Monitor Mount	FPMA103G	\$410

5033 Rear mounted Assistant utility	5033	\$1040
-------------------------------------	------	--------

**Optional Features:**

		<u>Part#</u>		
Touchpad for chair control	5000:	Chair mount	2000-230	\$426
		Delivery mount	2000-233	\$626
	2001:	Chair mount	2000-231	\$426
		Delivery mount	2000-234	\$626
Swivel instrument tray holder		4000-013-1	\$186	
Additional umbilical per foot		4000-054	\$32	
Additional handpiece station		4000-050	\$186	
Water/Air quick disconnect		4091-6	\$38	
Additional electronics cable pre-wiring		4000-180	\$76	
Wet/dry foot control (in lieu of)		4003-121	\$80	
Additional HVE, holder and tubing		4000-051	\$88	
Additional syringe, holder and tubing		4000-053	\$138	
Additional SE, holder and tubing		4000-052	\$88	
Telescoping assistant arm		4027	\$394	
Remote bowl flush and cup fill switch		4093-1	\$138	
Internal Umbilical <i>(for Westar chairs only)</i>		4000-054-1	\$146	
City/Bottle upgrade <i>(includes water MSO and routing valve)</i>		4095-7	\$204	

Standard Features:

4085 and 4081 Cart Systems:

- Three handpiece automatic valves
- Tri-Block valve
- Master on/off
- Individual water coolant adjustment
- Three way syringe with autoclavable tip
- Non-retracting water coolant
- Handpiece purge toggle
- Asepsis Midwest 4 hole tubing
- Standard foot control
- Utility Center
- Regulators, Filters, and Air/Water shut off valves
- Clean water system *(bottle water only)*
- Stainless Instrument Tray

Additional Features of the 4081 Dual Purpose Cart:

- Swing mount delivery for left/right convertibility
- HVE and SE
- Solids collector
- 1/4 water/air quick disconnect

Options for 4085 Doctor's Cart:

- 1/4" water/air quick disconnect 4091-6 \$ 38
- City/Bottle upgrade 4095-7 \$204

**4085 Doctor's Cart - \$3138**



**4089 Doctor's and Assistant Cart - \$3350**  
*(add vacuum instrumentation to 4085 Cart)*

**4081 Dual Purpose Cart - \$3716**



**4082 Assistant Cart - \$2440**  
*(Shown with optional Stainless Tray)*



Standard Features:

4082 Assistant Cart

- Swing mount arm for left/right convertibility
- HVE and SE
- Solids collector
- Clean water system *(bottle water only)*

Options for 4082 Assistant Cart:

- Additional HVE, hose and holder 4000-051 \$88
- 1/4" water/air quick disconnect 4091-6 \$38
- Stainless Instrument Tray 4003-102 \$60
- Plastic Instrument Tray 4003-237 \$16